



**BYD
CARE** 

The 'BYD CARE' text is in a bold, black, sans-serif font. To the right of the text is a red heart icon with a white outline, symbolizing care and health.

PRODUCT INTRODUCTION

**BYD N95 Respirator
(Surgical Mask)
GB 19038-2010**

N95 Respirator (Surgical Mask)



Foldable Non-sterile (LATEX)

BYD CARE



* Picture is for reference only, actual product may have differences.

- + Very close facial **fit**
- + \geq **95%** filtrations efficiency against solid and liquid aerosols free of oil.

Description

The N95 Respirator (Surgical Mask) is designed to help provide respiratory protection for the user. This product has a filter efficiency of 95% or greater against solid and liquid aerosols free of oil. This product contains no components made from natural rubber latex.

Use instructions

- Suitable for protecting the mouth, nose, and chin to airborne particles, and block liquid aerosols, blood, body fluids, and secretions
- For adult use only.
- Do not wash.
- Store in a cool, dry, clean place away from fire and contamination.
- Store between -20°C and 38°C (-4°F and 100°F) with relative humidity below 80%.
Expiration date: 2 years after production.

3-Step check before use

1. Do not use if the packaging has been opened or damaged.
2. Check if the elastic bands are in good condition and are not damaged.
3. Check that the metallic strip is not broken.

Wearing instructions

- 1 Unfold the mask and hold it with both hands, with the metallic strip facing up.
- 2 Place elastic bands around the neck and head respectively. Put the mask against your face covering both nose and mouth.
- 3 Adjust the metallic strip over the bridge of the nose by using two fingers to press down until achieving a close fit.
- 4 Perform a fit test according to instructions.



Fit test


1. Place both hands over the mask.
2. Take a deep breath and hold your breath for a few seconds, making sure the mask collapses inward as you inhale.
3. Then exhale and hold your breath for a few seconds, making sure the mask bulges outward as you exhale.
4. If air leaks, reposition the mask, tighten the metallic strip and reposition the elastic bands for a better fit. Follow steps 1-3 again until a tight seal has been achieved.
5. Masks that have passed the fitting test are now safe to use.
6. Facial hair or other facial features may prevent the mask from achieving a proper fit. If a tight seal is not achievable, please DO NOT use this product.
7. If you experience difficulty breathing, dizziness, or other symptoms, leave the contaminated area immediately.

N95 Respirator(Surgical Mask)



Foldable (Non-sterile) (Latex)

Specifications



Name	Parameter
Product name	N95 Respirator (Surgical Mask)
Material	Polypropylene spunbond nonwoven fabric (for surgical use), Polypropylene melt-blown nonwoven (for surgical use), hot air cotton, Polypropylene spunbond nonwoven fabric (for surgical use), Polyester and spandex fabric elastic bands, Aluminum nose clip
Product model	DE2322
Packaging specification	1pc per bag; 50pcs per box; 1000pcs per carton

GB 19083-2010 standard

WARNINGS

1. The mask does not eliminate the risk of contracting any disease or infection.
2. Improper use may lead to illness and even death.
3. Use this product immediately after the package is opened.
4. DO NOT sleep while wearing the mask.
5. Avoid hand contact within the interior part of the mask.
6. DO NOT use masks if they expired.
7. For one-time use only. Dispose the mask according to regulations.
8. DO NOT use if the package is damaged.
9. This N95 respirator is authorized for distribution under FDA Emergency Use Authorization only for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of N95 respirators under section 564 of the Act, 21 U.S.C. § 360bbb-3, unless the authorization is terminated or revoked sooner



Time use limitation

If the mask becomes damaged, soiled, or breathing becomes difficult, leave the contaminated area and replace and refit the new mask.

Important notice

To the extent permitted by law, BYD shall not be liable for any loss or damage including any loss of business, loss of profits, or for any indirect, special, incidental or consequential loss or damage arising from reliance upon any information herein provided by BYD. Nothing in this statement will be deemed to exclude or restrict BYD's liability for death or personal injury arising from its negligence.